

## Peptidyl-prolyl cis-trans isomerase A

Catalog No: #AP87060



Package Size: #AP87060-1 50ug #AP87060-2 100ug #AP87060-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Peptidyl-prolyl cis-trans isomerase A                              |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                                       |
| Species Reactivity | Bovine   |
| Other Names        | PPIase A, Cyclophilin A, Cyclosporin A-binding protein, Rotamase A |
| Accession No.      | Uniprot ID:P62935  |
| Uniprot            | P62935   |
| GeneID             | 281418;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)             |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PPIA

Note: This product is for in vitro research use only