

Peroxisome proliferator-activated receptor gamma coactivator-related protein 1

Catalog No: #AP87061

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87061-1 50ug #AP87061-2 100ug #AP87061-3 1mg

Description

Product Name	Peroxisome proliferator-activated receptor gamma coactivator-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PGC-1-related coactivator,PRC,Kiaa0595
Accession No.	Uniprot ID:Q6NZN1
Uniprot	Q6NZN1
GeneID	226169;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pprc1

Note: This product is for in vitro research use only