

Vitamin K-dependent protein S

Catalog No: #AP87063



Package Size: #AP87063-1 50ug #AP87063-2 100ug #AP87063-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitamin K-dependent protein S
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PROS
Accession No.	Uniprot ID:P07224
Uniprot	P07224
GeneID	282006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PROS1

Note: This product is for in vitro research use only