

Protein disulfide-isomerase

Catalog No: #AP87065



Package Size: #AP87065-1 50ug #AP87065-2 100ug #AP87065-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PDI, Cellular thyroid hormone-binding protein, Prolyl 4-hydroxylase subunit beta, p55, PDIA1
Accession No.	Uniprot ID:P05307
Uniprot	P05307
GeneID	281373;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: P4HB

Note: This product is for in vitro research use only