

Protein disulfide-isomerase A4

Catalog No: #AP87066



Package Size: #AP87066-1 50ug #AP87066-2 100ug #AP87066-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Endoplasmic reticulum resident protein 70,ER protein 70,Endoplasmic reticulum resident protein 72,ER protein 72,ERP70,ERP72
Accession No.	Uniprot ID:P13667
Uniprot	P13667
GeneID	9601;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDIA4

Note: This product is for in vitro research use only