

Peroxiredoxin-1

Catalog No: #AP87068



Package Size: #AP87068-1 50ug #AP87068-2 100ug #AP87068-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Natural killer cell-enhancing factor A,NKEF-A,Proliferation-associated gene protein,PAG,Thioredoxin peroxidase 2,Thioredoxin-dependent peroxide reductase 2,PAGA,PAGB,TDPX2
Accession No.	Uniprot ID:Q06830
Uniprot	Q06830
GeneID	5052;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRDX1

Note: This product is for in vitro research use only