

Prominin-1

Catalog No: #AP87070



Package Size: #AP87070-1 50ug #AP87070-2 100ug #AP87070-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prominin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Antigen AC133 homolog,Prominin-like protein 1,CD133,Prom,Prom1
Accession No.	Uniprot ID:O54990
Uniprot	O54990
GeneID	19126;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prom1

Note: This product is for in vitro research use only