

Parathyroid hormone-related protein

Catalog No: #AP87072



Package Size: #AP87072-1 50ug #AP87072-2 100ug #AP87072-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Parathyroid hormone-related protein |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | PTH-rP,Parathyroid hormone-like protein,PLP,Pthrp |
| Accession No. | Uniprot ID:P13085 |
| Uniprot | P13085 |
| GeneID | 24695; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pthlh

Note: This product is for in vitro research use only