

Pentraxin-related protein PTX3

Catalog No: #AP87073



Package Size: #AP87073-1 50ug #AP87073-2 100ug #AP87073-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pentraxin-related protein PTX3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pentraxin-related protein PTX3, Tumor necrosis factor-inducible gene 14 protein, TSG-14, Tsg14
Accession No.	Uniprot ID:P48759
Uniprot	P48759
GeneID	19288;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ptx3

Note: This product is for in vitro research use only