

# Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein

Catalog No: #AP87075

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87075-1 50ug #AP87075-2 100ug #AP87075-3 1mg

## Description

Product Name	Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 2 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	P-Rex2,DEP domain-containing protein 2,DEPDC2
Accession No.	Uniprot ID:Q70Z35
Uniprot	Q70Z35
GeneID	80243;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PREX2

Note: This product is for in vitro research use only