

# Peroxisome proliferator-activated receptor gamma coactivator 1-alpha

Catalog No: #AP87076

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87076-1 50ug #AP87076-2 100ug #AP87076-3 1mg

## Description

Product Name	Peroxisome proliferator-activated receptor gamma coactivator 1-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	PGC-1-alpha,PGC1,PGC1A,PPARGC1
Accession No.	Uniprot ID:Q865B6
Uniprot	Q865B6
GeneID	397013;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPARGC1A

Note: This product is for in vitro research use only