

DNA primase large subunit

Catalog No: #AP87077



Package Size: #AP87077-1 50ug #AP87077-2 100ug #AP87077-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA primase large subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	DNA primase 58 kDa subunit,p58
Accession No.	Uniprot ID:P33610
Uniprot	P33610
GeneID	19076;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prim2

Note: This product is for in vitro research use only