

Regenerating islet-derived protein 3-beta

Catalog No: #AP87080



Package Size: #AP87080-1 50ug #AP87080-2 100ug #AP87080-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Regenerating islet-derived protein 3-beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	REG-3-beta,Pancreatitis-associated protein 1,Regenerating islet-derived protein III-beta,Reg III-beta,Pap,Pap1
Accession No.	Uniprot ID:P35230
Uniprot	P35230
GeneID	18489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Reg3b

Note: This product is for in vitro research use only