

Polypyrimidine tract-binding protein 1

Catalog No: #AP87081



Package Size: #AP87081-1 50ug #AP87081-2 100ug #AP87081-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Polypyrimidine tract-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PTB,Heterogeneous nuclear ribonucleoprotein I,hnRNP I,Pyrimidine-binding protein,PYBP,Ptb,Pybp,Tbfii
Accession No.	Uniprot ID:Q00438
Uniprot	Q00438
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptbp1

Note: This product is for in vitro research use only