

Profilin-2

Catalog No: #AP87082



Package Size: #AP87082-1 50ug #AP87082-2 100ug #AP87082-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Profilin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Profilin II
Accession No.	Uniprot ID:Q9EPC6
Uniprot	Q9EPC6
GeneID	100909840;81531;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pfn2

Note: This product is for in vitro research use only