

Vitamin K-dependent protein C

Catalog No: #AP87086



Package Size: #AP87086-1 50ug #AP87086-2 100ug #AP87086-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vitamin K-dependent protein C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Anticoagulant protein C, Autoprothrombin IIA, Blood coagulation factor XIV
Accession No.	Uniprot ID:Q9GLP2
Uniprot	Q9GLP2
GeneID	396954;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PROC

Note: This product is for in vitro research use only