

Ras-related protein Rap-1A

Catalog No: #AP87089



Package Size: #AP87089-1 50ug #AP87089-2 100ug #AP87089-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras-related protein Rap-1A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ras-related protein Krev-1,Krev-1
Accession No.	Uniprot ID:P62835
Uniprot	P62835
GeneID	109905;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rap1a

Note: This product is for in vitro research use only