

Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C



Catalog No: #AP87090

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87090-1 50ug #AP87090-2 100ug #AP87090-3 1mg

Description

Product Name	Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cam-PDE 1C,hCam-3
Accession No.	Uniprot ID:Q14123
Uniprot	Q14123
GeneID	5137;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDE1C

Note: This product is for in vitro research use only