

## Podoplanin

Catalog No: #AP87092



Package Size: #AP87092-1 50ug #AP87092-2 100ug #AP87092-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Podoplanin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	E11 antigen epitope,Glycoprotein 38,Gp38,OTS-8,RTI140,T1-alpha,T1A,Type I cell 40 kDa protein,Gp38,Ots8
Accession No.	Uniprot ID:Q64294
Uniprot	Q64294
GeneID	54320;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pdpn

Note: This product is for in vitro research use only