

RNA-binding protein 3

Catalog No: #AP87093



Package Size: #AP87093-1 50ug #AP87093-2 100ug #AP87093-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RNA-binding protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RNA-binding motif protein 3
Accession No.	Uniprot ID:O89086
Uniprot	O89086
GeneID	19652;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rbm3

Note: This product is for in vitro research use only