

Prorelaxin

Catalog No: #AP87095



Package Size: #AP87095-1 50ug #AP87095-2 100ug #AP87095-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prorelaxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Accession No.	Uniprot ID:Q9TRM8
Uniprot	Q9TRM8
GeneID	403742;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RLN

Note: This product is for in vitro research use only