

Prohibitin-2

Catalog No: #AP87097



Package Size: #AP87097-1 50ug #AP87097-2 100ug #AP87097-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prohibitin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	B-cell receptor-associated protein BAP37,BAP-37,Bcap37
Accession No.	Uniprot ID:Q5XIH7
Uniprot	Q5XIH7
GeneID	114766;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Phb2

Note: This product is for in vitro research use only