

Response gene to complement 32 protein

Catalog No: #AP87099



Package Size: #AP87099-1 50ug #AP87099-2 100ug #AP87099-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Response gene to complement 32 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	RGC-32
Accession No.	Uniprot ID:Q9Z2P4
Uniprot	Q9Z2P4
GeneID	117183;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Rgc32

Note: This product is for in vitro research use only