

## L-seryl-tRNA(Sec) kinase

Catalog No: #AP87102



Package Size: #AP87102-1 50ug #AP87102-2 100ug #AP87102-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	L-seryl-tRNA(Sec) kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	O-phosphoseryl-tRNA(Sec) kinase
Accession No.	Uniprot ID:Q8BP74
Uniprot	Q8BP74
GeneID	214580;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Pstk

Note: This product is for in vitro research use only