

Progesterone receptor

Catalog No: #AP87103



Package Size: #AP87103-1 50ug #AP87103-2 100ug #AP87103-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Progesterone receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PR,Nuclear receptor subfamily 3 group C member 3,NR3C3
Accession No.	Uniprot ID:P06401
Uniprot	P06401
GeneID	5241;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PGR

Note: This product is for in vitro research use only