

Renin

Catalog No: #AP87107



Package Size: #AP87107-1 50ug #AP87107-2 100ug #AP87107-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Renin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Angiotensinogenase
Accession No.	Uniprot ID:Q6DYE7
Uniprot	Q6DYE7
GeneID	403838;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: REN

Note: This product is for in vitro research use only