

Junction plakoglobin

Catalog No: #AP87117



Package Size: #AP87117-1 50ug #AP87117-2 100ug #AP87117-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Junction plakoglobin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Catenin gamma,Desmoplakin III,Desmoplakin-3,CTNNG,DP3
Accession No.	Uniprot ID:P14923
Uniprot	P14923
GeneID	3728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:JUP

Note: This product is for in vitro research use only