

# 1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1

Catalog No: #AP87118

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP87118-1 50ug #AP87118-2 100ug #AP87118-3 1mg

## Description

Product Name	1-phosphatidylinositol 4,5-bisphosphate phosphodiesterase gamma-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PLC-148,Phosphoinositide phospholipase C-gamma-1,Phospholipase C-II,PLC-II,Phospholipase C-gamma-1,PLC-gamma-1
Accession No.	Uniprot ID:P08487
Uniprot	P08487
GeneID	281987;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PLCG1

Note: This product is for in vitro research use only