

Receptor-type tyrosine-protein phosphatase F

Catalog No: #AP87119



Package Size: #AP87119-1 50ug #AP87119-2 100ug #AP87119-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Receptor-type tyrosine-protein phosphatase F
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Leukocyte common antigen related,LAR,Lar
Accession No.	Uniprot ID:A2A8L5
Uniprot	A2A8L5
GeneID	19268;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ptpnf

Note: This product is for in vitro research use only