

# Prion-like protein doppel

Catalog No: #AP87131



Package Size: #AP87131-1 50ug #AP87131-2 100ug #AP87131-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prion-like protein doppel
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PrPLP,DPL
Accession No.	Uniprot ID:Q9GK16
Uniprot	Q9GK16
GeneID	281426;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PRND

Note: This product is for in vitro research use only