

Prostasin

Catalog No: #AP87132



Package Size: #AP87132-1 50ug #AP87132-2 100ug #AP87132-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prostasin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Channel-activating protease 1,CAP1,Serine protease 8
Accession No.	Uniprot ID:Q9ES87
Uniprot	Q9ES87
GeneID	192107;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prss8

Note: This product is for in vitro research use only