

## Ras-related protein Rab-8A

Catalog No: #AP87134



Package Size: #AP87134-1 50ug #AP87134-2 100ug #AP87134-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ras-related protein Rab-8A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Oncogene c-mel,Mel,Rab8
Accession No.	Uniprot ID:P55258
Uniprot	P55258
GeneID	17274;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Rab8a

Note: This product is for in vitro research use only