

Proteasome subunit beta type-8

Catalog No: #AP87137

Package Size: #AP87137-1 50ug #AP87137-2 100ug #AP87137-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteasome subunit beta type-8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Low molecular mass protein 7,Macropain subunit C13,Multicatalytic endopeptidase complex subunit C13,Proteasome component C13,Proteasome subunit beta-5i,Really interesting new gene 10 protein,LMP7,PSMB5i,RING10,Y2
Accession No.	Uniprot ID:P28062
Uniprot	P28062
GeneID	5696;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PSMB8

Note: This product is for in vitro research use only