

# Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform

Catalog No: #AP87143

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP87143-1 50ug #AP87143-2 100ug #AP87143-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phosphatidylinositol 4,5-bisphosphate 3-kinase catalytic subunit delta isoform
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PI3-kinase subunit delta, Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit delta, PtdIns-3-kinase subunit p110-delta
Accession No.	Uniprot ID:O35904
Uniprot	O35904
GeneID	18707;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Pik3cd

Note: This product is for in vitro research use only