

Prostaglandin G/H synthase 1

Catalog No: #AP87144



Package Size: #AP87144-1 50ug #AP87144-2 100ug #AP87144-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin G/H synthase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cyclooxygenase-1,COX-1,Prostaglandin H2 synthase 1,PGH synthase 1,Prostaglandin-endoperoxide synthase 1,COX1
Accession No.	Uniprot ID:P23219
Uniprot	P23219
GeneID	5742;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTGS1

Note: This product is for in vitro research use only