

Retinoic acid receptor alpha

Catalog No: #AP87145



Package Size: #AP87145-1 50ug #AP87145-2 100ug #AP87145-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RAR-alpha, Nuclear receptor subfamily 1 group B member 1, NR1B1
Accession No.	Uniprot ID:Q90966
Uniprot	Q90966
GeneID	395213;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RARA

Note: This product is for in vitro research use only