

Retinoic acid receptor beta

Catalog No: #AP87146



Package Size: #AP87146-1 50ug #AP87146-2 100ug #AP87146-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinoic acid receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	RAR-beta,HBV-activated protein,Nuclear receptor subfamily 1 group B member 2,RAR-epsilon,HAP,NR1B2
Accession No.	Uniprot ID:P10826
Uniprot	P10826
GeneID	5915;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:RARB

Note: This product is for in vitro research use only