

Basement membrane-specific heparan sulfate proteoglycan core protein

Catalog No: #AP87147

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP87147-1 50ug #AP87147-2 100ug #AP87147-3 1mg

Description

Product Name	Basement membrane-specific heparan sulfate proteoglycan core protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HSPG
Accession No.	Uniprot ID:Q05793
Uniprot	Q05793
GeneID	15530;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hspg2

Note: This product is for in vitro research use only