

Protein S100-A4

Catalog No: #AP87148



Package Size: #AP87148-1 50ug #AP87148-2 100ug #AP87148-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Metastasin,S100 calcium-binding protein A4,MTS1
Accession No.	Uniprot ID:Q9TV56
Uniprot	Q9TV56
GeneID	403787;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100A4

Note: This product is for in vitro research use only