

## Sphingosine 1-phosphate receptor 2

Catalog No: #AP87149



Package Size: #AP87149-1 50ug #AP87149-2 100ug #AP87149-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sphingosine 1-phosphate receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	S1P receptor 2, Endothelial differentiation G-protein coupled receptor 5, Lysophospholipid receptor B2, Sphingosine 1-phosphate receptor Edg-5, S1P receptor Edg-5, Edg5, Gpcr13, Lpb2
Accession No.	Uniprot ID:P52592
Uniprot	P52592
GeneID	14739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: S1pr2

Note: This product is for in vitro research use only