

Ribonucleoside-diphosphate reductase subunit M2

Catalog No: #AP87152



Package Size: #AP87152-1 50ug #AP87152-2 100ug #AP87152-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribonucleoside-diphosphate reductase subunit M2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Ribonucleotide reductase small chain,Ribonucleotide reductase small subunit
Accession No.	Uniprot ID:P11157
Uniprot	P11157
GeneID	20135;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rrm2

Note: This product is for in vitro research use only