

Phospholipase D2

Catalog No: #AP87154



Package Size: #AP87154-1 50ug #AP87154-2 100ug #AP87154-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase D2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PLD 2,Choline phosphatase 2,Phosphatidylcholine-hydrolyzing phospholipase D2
Accession No.	Uniprot ID:Q0V8L6
Uniprot	Q0V8L6
GeneID	522159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PLD2

Note: This product is for in vitro research use only