

Retinol-binding protein 1

Catalog No: #AP87158



Package Size: #AP87158-1 50ug #AP87158-2 100ug #AP87158-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retinol-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cellular retinol-binding protein,CRBP,Cellular retinol-binding protein I,CRBP-I,Crbpi,Rbp-1
Accession No.	Uniprot ID:Q00915
Uniprot	Q00915
GeneID	19659;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Rbp1

Note: This product is for in vitro research use only