

Prostaglandin F2-alpha receptor

Catalog No: #AP87160



Package Size: #AP87160-1 50ug #AP87160-2 100ug #AP87160-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prostaglandin F2-alpha receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PGF receptor, Prostanoid FP receptor
Accession No.	Uniprot ID:P37289
Uniprot	P37289
GeneID	282020;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PTGFR

Note: This product is for in vitro research use only