

Rho-related GTP-binding protein RhoC

Catalog No: #AP87163



Package Size: #AP87163-1 50ug #AP87163-2 100ug #AP87163-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rho-related GTP-binding protein RhoC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Rho cDNA clone 9,h9,ARH9,ARHC
Accession No.	Uniprot ID:P08134
Uniprot	P08134
GeneID	389;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: RHOC

Note: This product is for in vitro research use only