

Parathyroid hormone

Catalog No: #AP87164



Package Size: #AP87164-1 50ug #AP87164-2 100ug #AP87164-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Parathyroid hormone
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PTH,Parathormone,Parathyrin
Accession No.	Uniprot ID:P01270
Uniprot	P01270
GeneID	5741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTH

Note: This product is for in vitro research use only