

Pleiotrophin

Catalog No: #AP87169



Package Size: #AP87169-1 50ug #AP87169-2 100ug #AP87169-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pleiotrophin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PTN,Heparin-binding brain mitogen,HBBM,Heparin-binding growth factor 8,HBGF-8,Heparin-binding growth-associated molecule,HB-GAM,Heparin-binding neurite-promoting factor,p18
Accession No.	Uniprot ID:P21782
Uniprot	P21782
GeneID	280904;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PTN

Note: This product is for in vitro research use only