

2'-5'-oligoadenylate synthase 2

Catalog No: #AP87170

Package Size: #AP87170-1 50ug #AP87170-2 100ug #AP87170-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	2'-5'-oligoadenylate synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	(2-5')oligo(A) synthase 2,p69 OAS / p71 OAS,p69OAS / p71OAS
Accession No.	Uniprot ID:P29728
Uniprot	P29728
GeneID	4939;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OAS2

Note: This product is for in vitro research use only