

Peptide YY

Catalog No: #AP87172



Package Size: #AP87172-1 50ug #AP87172-2 100ug #AP87172-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peptide YY
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PYY,Peptide tyrosine tyrosine
Accession No.	Uniprot ID:P10631
Uniprot	P10631
GeneID	287730;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pyy

Note: This product is for in vitro research use only