

Prefoldin subunit 2

Catalog No: #AP87174



Package Size: #AP87174-1 50ug #AP87174-2 100ug #AP87174-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prefoldin subunit 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:B0BN18
Uniprot	B0BN18
GeneID	685607;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pfdn2

Note: This product is for in vitro research use only